

QUICK GUIDE to USING the VARIAN 300 NMR SPECTROMETERS (vnmr 6.1c)

To change sample & set ²H lock: **Acqi** button opens ACQUISITION menu. **Eject** the CDCl₃ standard. Carefully adjust tube position in sample holder (spinner) with depth gauge. Center short samples (< 5 cm) relative to the detection coil. **Insert** sample in magnet. Open the **LOCK** menu on the **Acqi** panel. Turn on SPIN (20 Hz; no spin for J-Young tubes). Turn **Lock off** to while adjusting **Z0** to match the deuterium resonance of the solvent.

CDCl₃ **INOVA-300:** **Z0= +39, lockpower=30, lockgain=45 (max 48)** lock level ~68%
MERCURY-300: **Z0=+387, lockpower=25, lockgain=26 (max 38)** lock level ~68%
GEMINI-300: **Z0=+232, lockpower=30, lockgain=20 (max 30)** lock level ~53%

These parameters differ on each spectrometer and for each deuterated solvent. Note that **Z0** drifts over a period of several months; When adjusting **Z0** it is helpful to temporarily increase the **lockpower** (be sure to turn it back down!). Engage the lock by clicking **Lock on**. *Autolock is not recommended.*

On the 300-NMR systems **lockphase** is only slightly sample dependent, but generally can be left alone. If necessary, tweak **lockphase** by +/-4 units at a time to optimize lock level. It may change after a power failure.

To shim: **Acqi** => **SHIM** (ACQUISITION menu)

Choose **manual** shimming. Open the **shim** menu with **RIGHT** mouse button and select "*fine Z*".

Alternately adjust **Z1** and **Z2** to maximize the lock level indicated by the color bars in the shim window.

Start with increments of +/-64, then +/-16. For more critical shimming use increments of +/-4 and also adjust **Z3**.

Repeat the cycle of adjusting **Z1/Z2/Z3** to achieve a maximum lock level, usually about 55-65% for CDCl₃.

Exit the ACQUISITION menu by clicking the **CLOSE** button.

Note: shim and lock settings are shared from one user to the next.

If necessary, reload a fresh copy of the standard shim values by typing: **rts('std5') su**

To choose nucleus/solvent and modify/load parameters: **Main Menu** => **Setup** => **H1, CDCl3**

To list parameters, type **dg**, **dg1**, or **dgs** for *Acquisition Display*, or *Shim* parameter groups, respectively.

Change parameters values by typing on the command line, e.g.: **nt=16** to collect 16 scans, or **sw=10000** to set the spectral width to 10,000 Hz. To check a parameter value, type the parameter name followed by a "?", e.g., **sw?**

To acquire & process spectrum:

Type **ga** (*get acquisition*) in the command line to start acquisition. Upon completion, spectrum is processed and displayed automatically. Type **dscale** to view the ppm scale. Auto-phase spectrum with **aph0** or phase interactively with the **Phase** menu button. To change the line-broadening and re-process, type: **lb=1 wft**

To expand the spectrum **Main Menu** => **Display** => **Interactive**

Set cursors on left and right of the spectrum with left/right mouse buttons and click on **Expand**.

To scale to the biggest peak type **vsadj**. Change the vertical expansion by dragging a peak up/down with the middle mouse button or with **vs** parameter; e.g., **vs=vs*10** increases the vertical expansion by tenfold.

To show the full spectrum, type **f full** or click the **Full** button. **Resize** increases/decreases the window size.

To display partial (segmented) integrals: **Main Menu** => **Display** => **Interactive** => **Part Integral**

Type the command **cz** to clear old partial integral marks. Click **Lvl/tlt** to adjust baseline tilt/offset.

Click on **Resets** and set new marks w/left mouse button. Use **isadj** to scale and **Set Int** to normalize integrals.

To pick peaks: **Main Menu** => **Display** => **Interactive** => **Th**

Set vertical threshold to desired level (yellow horizontal line) with mouse.

Type **dpf** in command line to view chemical shifts of the picked peaks.

To print:

For a basic plot, enter: **pl pscale page**. Plots are built up by a series of commands, terminated by **page**.

Optional commands: **pap** (*all parameters*), **ppa** (*parameters in "English"*), **pir** (*integrals*) and **ppf** (*picked peaks*).

To save data as a disk file:

To store a spectrum (including parameters, shims, text, etc) on the file server, type **svf('vnmrdata/filename')**.

To retrieve a spectrum file, type **rt('vnmrdata/filename')** or use the **File** => **Load** menu buttons.